


Angled Screw Abutment

Additional Option for a Successful Dentium Implant Restoration



Simple & Predictable Over 15 Years of Clinical Evidence

over 15 years
of Long term data



06. 20. 2002
Pre-op

07. 02. 2002
Post-op

11. 29. 2002
Final Prosthesis

06. 02. 2004
1 Years 2 months

11. 15. 2007
5 Years

07. 24. 2009
6 Years 8 months

05. 23. 2012
9 Years 6 months

11. 25. 2017
15 Years




DentiumUSA
Developed by Clinicians for Clinicians

6731 Katella Ave. Cypress, CA 90630
Tel. 1-877-304-6752 Fax. 1-714-226-0019
cs@dentiumusa.com

PN2997-2

Angled Screw Abutment



Screw Abutment



REF	Diameter	Angle	G/H (mm)		Type
			Short	Long	
ASA45151018HT	Ø4.5	15°	1.0	2.2	Hex
ASA45151018N			1.0	2.2	Non-Hex
ASA45152018HT			2.0	3.2	Hex
ASA45152018N			2.0	3.2	Non-Hex
ASA45301018HT		30°	1.0	3.3	Hex
ASA45301018N			1.0	3.3	Non-Hex
ASA45302018HT			2.0	4.3	Hex
ASA45302018N			2.0	4.3	Non-Hex
ASA55151518HT	Ø5.5	15°	1.5	3.0	Hex
ASA55151518N			1.5	3.0	Non-Hex
ASA55153018HT			3.0	4.5	Hex
ASA55153018N			3.0	4.5	Non-Hex
ASA55301518HT		30°	1.5	4.3	Hex
ASA55301518N			1.5	4.3	Non-Hex
ASA55303018HT			3.0	5.5	Hex
ASA55303018N			3.0	5.5	Non-Hex

REF	Diameter	G/H(mm)
SAB4510L	Ø4.5	1.0
SAB4515L		1.5
SAB4525L		2.5
SAB4535L		3.5
SAB4545L		4.5
SAB4555L		5.5
SAB5515L	Ø5.5	1.5
SAB5525L		2.5
SAB5535L		3.5
SAB5545L		4.5
SAB5555L		5.5

Related Components



Cylinder

Impression Coping

Lab Analog



DentiumUSA
Developed by Clinicians for Clinicians

6731 Katella Ave. Cypress, CA 90630
Tel. 1-877-304-6752 Fax. 1-714-226-0019
cs@dentiumusa.com